

REXMONO[®]

Synthetic Absorbable
POLYGLECAPRONE
SUTURES






















REXMONO[®] SUTURES

REXMONO[®] suture is fast absorbable synthetic monofilament suture composed of a Poly(glycolide-co-caprolactone) and is available both dyed and undyed. When dyed, only FDA-approved color additives such as D&C Violet No. 2 are used. The undyed suture is beige in color. REXMONO meets all the requirements of the European Pharmacopoeia (EP) for Synthetic Absorbable Monofilament Sutures. The suture is smooth, making it easier to pass through the tissues as compared to braided suture. The suture is rarely used for percutaneous skin closure and is not used in areas of high tension. REXMONO has low tissue reactivity and maintains 50 to 60% in-vivo tensile strength at 1 weeks. Complete absorption in tissues takes around 90 to 110 days.














Needle Point Type

- Taper point ⦿
- Reverse Cutting ▼
- Conventional cutting ▲
- Spatula ▽
- Blunt ⊙
- Taper Cutting ⊗

Needle Type	Length & Color	Suture U.S.P (Metric)				
		6-0 (0.7)	5-0 (1)	4-0 (1.5)	3-0 (2)	2-0 (3)
8mm						
 1/4	 50cm Violet	MA608FD		MA408FD		
	3/8	MB608D				
 3/8	 50cm Violet	MM608				
9mm						
 3/8	 50cm Violet	MB609				
 1/2	 30cm Violet	MN609	MN509			
 3/8	 50cm Violet		MB509			
11mm						
 3/8	 50cm Violet	ME6111	ME511	ME411		
	 50cm Violet	ME611D				
	1/2	MF611	MF511			
 1/2	 30cm Violet	MN611	MN511	MN411		
	str.			MS411D		
13mm						
 3/8	 45cm Violet		ME513S			
	50cm Violet	ME613	ME513	ME413		
	1/2	MF613	MF513			
 1/2	 70cm Violet	MN613	MN513			
	3/8			MM413		

REXMONO[®]








POLYGLYCAPRONE (Glycolide-co-caprolactone) - Monofilament










Needle Type	Length & Color	Suture U.S.P (Metric)				
		6-0 (0.7)	5-0 (1)	4-0 (1.5)	3-0 (2)	2-0 (3)
14mm						
⊙ 1/2 	30cm Violet	MN614	MN514			
	50cm Violet			MN414		
15mm						
▽ 3/8 	50cm Violet	ME615	ME515	ME415	ME315	
	50cm Violet 	ME615D		ME415D		
1/2 	50cm Violet		MF515			
⊙ 3/8 	75cm Violet			MM415		
17mm						
⊙ 1/2 	50cm Violet		MN517S	MN417S	MN317S	MN217
	75cm Violet	MN617	MN517	MN417	MN317	
⊕ 5/8 	50cm Violet			MT417		
▽ 3/8 	50cm Violet		ME517			
	75cm Violet			ME417	ME317	
1/2 	50cm Violet		MF517	MF417		
18mm						
⊙ 1/2 	50cm Violet	MN618	MN518	MN418	MN318	
▽ 3/8 	50cm Violet			ME418	ME318	
1/2 	50cm Violet		MF518	MF418	MF318	
19mm						
▽ 3/8 	50cm Violet		ME519	ME419	ME319	

Needle Type	Length & Color	Suture U.S.P (Metric)				
		5-0 (1)	4-0 (1.5)	3-0 (2)	2-0 (3)	0 (3.5)
19mm						
⊙ 3/8 	50cm Violet	MM519	MM419	MM319		
1/2 	50cm Violet	MN519	MN419	MN319		
21mm						
▽ 3/8 	50cm Violet		ME421	ME321	ME221	
1/2 	50cm Violet	MF521	MF421	MF321	MF221	
⊙ 1/2 	50cm Violet 75cm Violet	MN521				
			MN421	MN321	MN221	MN021
22mm						
⊙ 1/2 	75cm Violet		MN422	MN322	MN222	
23mm						
⊙ 1/2 	75cm Violet		MN423	MN323	MN223	
24mm						
⊙ 1/2 	75cm Violet 50cm Violet		MN424		MN324	
25mm						
▽ 3/8 	75cm Violet		ME425	ME325	ME225	
⊙ 1/2 	75cm Violet		MN425	MN325	MN225	MN025
5/8 	75cm Violet			MN325F		
⊕ 1/2 	75cm Violet				MT225	

REXMONO[®]


POLYGLCAPRONE (Glycolide-co-caprolactone) - Monofilament



Needle Type	Length & Color	Suture U.S.P (Metric)				
		4-0 (1.5)	3-0 (2)	2-0 (3)	0 (3.5)	1 (4)
26mm						
⊙ 5/8 	75cm Violet		MN326F	MN226F		
⊙ 1/2 	75cm Violet	MN426	MN326S	MN226S	MN026	
REMOVABLE 	90cm Violet			MN226		
	45cm Violet		MN326C (8strands)			
▽ 3/8 	75cm Violet		ME326	ME226	ME026	
⊕ 1/2 	75cm Violet			MT226		
28mm						
⊙ 1/2 	75cm Violet	MN428	MN328	MN228	MN028	
29mm						
▽ 3/8 	75cm Violet		ME329			
1/2 	75cm Violet		ME329	ME229		
⊙ 3/8 	75cm Violet			MM229	MM029	
1/2 	75cm Violet		MN329	MN229	MN029	MN0129

Needle Type	Length & Color	Suture U.S.P (Metric)				
		3-0 (2)	2-0 (3)	0 (3.5)	1 (4)	2 (5)
29mm						
⊙ 1/2 	75cm Violet		MN229D			
⊕ 1/2 	75cm Violet		MT229		MT0129	
30mm						
▽ 3/8 	75cm Violet			ME030		
⊙ 1/2 	75cm Violet			MN030		
31mm						
⊙ 1/2 	75cm Violet	MN331	MN231	MN031	MN0131	
33mm						
⊙ 1/2 	75cm Violet	MN333	MN233	MN033	MN0133	
36mm						
⊙ 1/2 	75cm Violet		MN236	MN036	MN0136	
	90cm Violet			MN036L		
5/8 	75cm Violet		MA236F			
⊕ 1/2 	75cm Violet		MT236	MT036		

REXMONO®







POLYGLCAPRONE (Glycolide-co-caprolactone) - Monofilament






Needle Type	Length & Color	Suture U.S.P (Metric)				
		3-0 (2)	2-0 (3)	0 (3.5)	1 (4)	2 (5)
36mm						
⊕ 1/2 	75cm Violet				MT0136	MT0236
▽ 3/8 	75cm Violet				ME0136	
37mm						
⊕ 1/2 	75cm Violet	MN337	MN237	MN037	MN0137	
⊕ 5/8 	75cm Violet			MN037F		
▽ 3/8 	75cm Violet		ME237			
⊕ 1/2 	75cm Violet		MT237	MT037	MT0137	
38mm						
⊕ 1/2 	75cm Violet	MN338	MN238	MN038		
40mm						
⊕ 1/2 	90cm Violet				MN0140	
41mm						
⊕ 3/8 	75cm Violet		MM241	MM041		

Needle Type	Length & Color	Suture U.S.P (Metric)				
		3-0 (2)	2-0 (3)	0 (3.5)	1 (4)	2 (5)
41mm						
⊙ 1/2 	75cm Violet		MN241	MN041	MN0141	MN0241
▽ 3/8 	75cm Violet			ME041		
1/2	75cm Violet			MF041	MF0141	MF0241
⊙ 1/2	60cm Violet					MT0141D
	75cm Violet					MT0141
⊙ 3/8	75cm Violet					MU0141
43mm						
⊙ 1/2	75cm Violet		MN243	MN043	MN0143	
45mm						
⊙ 3/8	75cm Violet		MM205	MM045	MM0145	

REXMONO[®]





POLYGLYCAPRONE (Glycolide-co-caprolactone) - Monofilament


Needle Type	Length & Color	Suture U.S.P (Metric)			
		2-0 (3)	0 (3.5)	1 (4)	2 (5)
45mm					
 1/2 	75cm Violet	MN245	MN045	MN0145	MN0245
46mm					
 1/2 	75cm Violet	MT246		MT0146	
 1/2 	75cm Violet		MN046	MN0146	MN0246
 1/2 	75cm Violet		MF046	MF0146	
48mm					
 1/2 	90cm Violet	MN248		MN0148	

Needle Type	Length & Color	Suture U.S.P (Metric)			
		2-0 (3)	0 (3.5)	1 (4)	2 (5)
48mm					
⊙ 1/2 	75cm Violet	MN248S	MN048	MN0148S	
▽ 1/2 	75cm Violet				MF0248
49mm					
⊙ 1/2 	75cm Violet	MN249	MN049	MN0149	MN0249
50mm					
⊙ str. 	75cm Violet	MS250			
51mm					
▽ 1/2 	75cm Violet		MF051	MF0151	

REXMONO[®]

POLYGLYCAPRONE (Glycolide-co-caprolactone) - Monofilament









Needle Type	Length & Color	Suture U.S.P (Metric)			
		2-0 (3)	0 (3.5)	1 (4)	2 (5)
56mm 	75cm Violet	MN256	MN056		
63mm 	75cm Violet		MU063	MU0163	
66mm 	75cm Violet	MM266	MM066		
	75cm Violet		MN066	MN0166	MN0266

Needle Type	Length & Color	Suture U.S.P (Metric)			
		2-0 (3)	0 (3.5)	1 (4)	2 (5)
85mm <input checked="" type="radio"/> 3/8 	75cm Violet			MU0185	

REXMONO[®]

POLYGLCAPRONE (Glycolide-co-caprolactone) - Monofilament

Needle Type	Length & Color	Suture U.S.P (Metric)				
		6-0 (0.7)	5-0 (1)	4-0 (1.5)	3-0 (2)	2-0 (3)
8mm						
 1/4 	45cm Undyed	MWA608F				
11mm						
 3/8 	45cm Undyed	MWE611	MWE511			
13mm						
 3/8 	45cm Undyed	MWE613	MWE513	MWE413		
16mm						
 3/8 	45cm Undyed		MWE516	MWE416		
17mm						
 1/2 	75cm Undyed				MWN317	
19mm						
 3/8 	45cm Undyed		MWE519	MWE419S	MWE319S	
	75cm Undyed			MWE419	MWE319	

Needle Type	Length & Color	Suture U.S.P (Metric)				
		6-0 (0.7)	5-0 (1)	4-0 (1.5)	3-0 (2)	2-0 (3)
24mm						
 3/8 	45cm Undyed			MWE424		
	75cm Undyed			MWE424S	MWE324	MWE224
26mm						
 3/8 	75cm Undyed			MWE426	MWE326	
48mm						
 str. 	75cm Undyed			MWS448		
60mm						
 str. 	75cm Undyed			MWS460	MWS360	



(주)에스엠이엔지
SM Eng
www.sm-eng.net

46, Nakdong-daero 1302beon-gil, Sasang-gu, Busan, Korea
TEL : 82-51-305-8016 FAX : 82-51-305-8021
E-mail : sales@sm-eng.net <http://www.sm-eng.net>